





The Sanibel Supply Catalogue

Quality Supplies for your clinic



Contents



Page 5 - 7

Disposable ear tips



Page 9 - 10

Probe tips



Page 11 - 12

Supplies for insert earphones



Page 13

Sanibel™ Infant EarCup™



Page 14

REM probe tubes for HA fitting



Page 15 - 19

Disposable electrodes

List of abbreviations

- ABR** Auditory Brainstem Response
- AUD** Audiometry
- BAL** Balance
- HIF** Hearing Instrument Fitting
- IMP** Impedance
- OAE** Otoacoustic Emissions



Page 20 - 21

Electrode cables



Page 21-22

Patient preparation



Page 23

Protection & cleaning



Page 24 - 27

Parts for audiometry



Page 28

Balance supplies



Page 29

Paper



Page 30 - 31

Miscellaneous



Single patient use only!



Product items are not made with natural rubber latex

Why use ear tips for single use only?

Sanibel Supply® recommends using disposable ear tips for single use only to stop ear infections from spreading

At Sanibel Supply® we think that most audiologists would agree with us that it is crucial to avoid inadvertent transmission of disease within the audiology clinic. The good news is that it is quite easy to ensure that ear infections are not spread from one patient to another; use a new ear tip for each patient, and you will significantly reduce or eliminate the incidences of cross contamination.

Cerumen may be an infectious agent

Cerumen is not considered an infectious agent unless contaminated with blood or mucus, and due to its color and viscosity, visual detection of blood or ear discharge contaminants may be difficult. Therefore, you should treat cerumen as a potentially infectious agent. It is especially important to be mindful of using new tips when testing children and those with weaker immune systems, such as the frail and elderly as they are more susceptible to infections.

Are disposable ear tips expensive?

Sanibel Supply® offers volume discounts on ear tips. This will often offset the time and cost involved in the ear tip cleaning process.

Be more cost effective in your clinic

By using a new ear tip pr. every patient you reach a seal much more reliable than with re-used ear tips. How often have you not been fumbling for minutes before you get a seal with a patient? Using the right ear tip and a fresh one pr. patient ensures you a reliable fit with every patient, meaning you can calculate less minutes for testing a patient and in most cases allow you to have an extra patient pr. day without having stressed out personnel. The cost of using one ear tip pr. patient is far less than what you earn on your next patient, plus you can promote your clinic as taking patient safety and hygiene seriously.

6 strong reasons why you should use a new ear tip after each patient:

1

You reduce the risk of cross contamination between your patients

Unfortunately, it is not possible to be 100% certain that complete decontamination has taken place when cleaning ear tips.

2

You minimize the risk of non-sealing ear tips

Ear tips degrade over time, which can make it more difficult to attain a good seal, hence leading to less accurate results.

3

Your patients avoid allergic reactions from cleaning agent residue

Residue from cleaning agents may cause allergic reactions in some patients.

4

Your patients will appreciate it

Patients are more comfortable knowing that new ear tips are being used.

5

It's so much easier

Cleaning ear tips takes time, so why not give your staff the time to do what really matters?

6

It's so much faster

A fresh ear tip creates a seal quicker and more reliable, so you do not waste time performing a standard procedure.



Selecting the correct ear tip

Age	Ear tips inserted against ear canal	Ear tips inserted into ear canal	Ear tips for Echo-Screen®	Ear tips for AccuScreen®
0.5 y/o		 Green pediatric ear tip	 Ear tip 4.0 mm, ANA	 Ear tip 3.7 mm, AOT
		 Ear tip 3-5 mm, flanged	 Ear tip 4.5 mm, ANA	 Ear tip 4.0 mm, AOT
	 Ear tip 7.5 mm, umbrella	 Ear tip 4-7 mm, flanged	 Ear tip 5.0 mm, ANA	 Ear tip 4.5 mm, AOT
	 Ear tip 9 mm, umbrella	 Ear tip 5-8 mm, flanged	 Ear tip infant tree, ANA	 Ear tip 5.0 mm, AOT
	 Ear tip 10 mm, umbrella	 Ear tip 7 mm, mushroom	 Ear tip mushroom, ANA	 Ear tip 4-7 mm, AOT
	 Ear tip 11 mm, umbrella	 Ear tip 8 mm, mushroom		
10 y/o	 Ear tip 13 mm, umbrella	 Ear tip 9 mm, mushroom		
		 Ear tip 10 mm, mushroom		
Adults	 Ear tip 15 mm, umbrella	 Ear tip 11 mm, mushroom		
		 Ear tip 12 mm, mushroom		
	 Ear tip 19 mm, umbrella	 Ear tip 13 mm, mushroom		
		 Ear tip 14 mm, mushroom		
	 Ear tip 22 mm, umbrella	 Ear tip 15 mm, mushroom		
	 Ear tip 19 mm, mushroom			

Selecting the correct size of ear tip takes practice

In general, flanged or mushroom-shaped ear tips will be used for most testing (Impedance, DPOAE, TEOAE and Automated ABR).

Sometimes, an umbrella-shaped ear tip will be used for performing a quick impedance screening.

The table can be used as a guide to help you select the appropriate ear tip size.

Disposable ear tips



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags

ADI ear tips

Sanibel™ flanged infant ear tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip: Titan™, AT235™, AA222™, Eclipse™, Otowave™, easyTymp™, touchTymp™, Grason-Stadler® TympStar™ and new SPL60 (RECD).

Part no.	Description	Unit
8502838	Pediatric (green)	Bag/100 pcs
8012963	3-5 mm (red)	Bag/25 pcs
8012964	3-5 mm (red)	Bag/100 pcs
8012965	4-7 mm (blue)	Bag/25 pcs
8012966	4-7 mm (blue)	Bag/100 pcs
8012967	5-8 mm (yellow)	Bag/25 pcs
8012968	5-8 mm (yellow)	Bag/100 pcs



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags



ADI ear tips

Sanibel™ mushroom ear tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip: Titan™, AT235™, AA222™, Eclipse™, Otowave™, easyTymp™, touchTymp™, Grason-Stadler® TympStar™ and new SPL60 (RECD).

Part no.	Description	Unit
8013001	7 mm (blue)	Bag/25 pcs
8013002	7 mm (blue)	Bag/100 pcs
8013003	8 mm (yellow)	Bag/25 pcs
8013004	8 mm (yellow)	Bag/100 pcs
8012969	9 mm (green)	Bag/25 pcs
8012970	9 mm (green)	Bag/100 pcs
8012971	10 mm (red)	Bag/25 pcs
8012972	10 mm (red)	Bag/100 pcs
8012973	11 mm (blue)	Bag/25 pcs
8012974	11 mm (blue)	Bag/100 pcs
8012975	12 mm (yellow)	Bag/25 pcs
8012976	12 mm (yellow)	Bag/100 pcs
8012977	13 mm (green)	Bag/25 pcs
8012978	13 mm (green)	Bag/100 pcs
8012979	14 mm (red)	Bag/25 pcs
8012980	14 mm (red)	Bag/100 pcs
8012981	15 mm (blue)	Bag/25 pcs
8012982	15 mm (blue)	Bag/100 pcs
8012983	19 mm (yellow)	Bag/25 pcs
8012984	19 mm (yellow)	Bag/100 pcs



Medical grade silicone rubber
Biocompatible

VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags



ADI ear tips

Sanibel™ umbrella ear tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip: Titan™, AT235™, AA222™, Eclipse™, Otowave™, easyTymp™, touchTymp™, Grason-Stadler® TympStar™ and new SPL60 (RECD).

Part no.	Description	Unit
8012985	7.5 mm (green)	Bag/25 pcs
8012986	7.5 mm (green)	Bag/100 pcs
8012987	9 mm (red)	Bag/25 pcs
8012988	9 mm (red)	Bag/100 pcs
8012989	10 mm (blue)	Bag/25 pcs
8012990	10 mm (blue)	Bag/100 pcs
8012991	11 mm (yellow)	Bag/25 pcs
8012992	11 mm (yellow)	Bag/100 pcs
8012993	13 mm (green)	Bag/25 pcs
8012994	13 mm (green)	Bag/100 pcs
8012995	15 mm (red)	Bag/25 pcs
8012996	15 mm (red)	Bag/100 pcs
8012997	19 mm (blue)	Bag/25 pcs
8012998	19 mm (blue)	Bag/100 pcs
8012999	22 mm (yellow)	Bag/25 pcs
8013000	22 mm (yellow)	Bag/100 pcs



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags



Single use ear tips for immittance instruments & OAE screeners

Sanibel™ infant ear tips

For use with many acoustic immittance instruments and otoacoustic emission systems, including the MAICO® ERO-SCAN™, Interacoustics® OtoRead™ and Grason-Stadler® TympStar™.

Part no.	Description	Unit
8013005	3 mm (yellow)	Bag/25 pcs
8013006	3 mm (yellow)	Bag/100 pcs
8002036	4 mm (red)	Bag/25 pcs
8002035	4 mm flanged (red)	Bag/100 pcs





AKR Single use ear tips for immittance instruments & OAE screeners

Mushroom ear tips

For use with many acoustic immittance instruments and otoacoustic emission systems, including the MAICO® ERO•SCAN™, Interacoustics® OtoRead™ and Grason-Stadler® TympStar™.

Part no.	Description	Unit
8500560	6 mm (green)	Bag/100 pcs
8500550	7 mm (yellow)	Bag/100 pcs
8500540	8 mm (red)	Bag/100 pcs
8500530	9 mm (blue)	Bag/100 pcs
8500520	10 mm (green)	Bag/100 pcs
8500510	11 mm (yellow)	Bag/100 pcs
8500500	12 mm (red)	Bag/100 pcs
8500490	13 mm (blue)	Bag/100 pcs



Single use ear tips for immittance instruments & OAE screeners

Mushroom ear tips

For use with Grason-Stadler® TympStar Pro™ and Corti™, Interacoustics® AT235™, AA222™ and Titan™.

Part no.	Description	Unit
8109426	3-18.5 mm	Sample kit/25 pcs
8109427	3-18.5 mm	Sample kit/10 pcs
8109428	3-6 mm (green)	Bag/100 pcs
8109429	6 mm (yellow)	Bag/100 pcs
8109430	7 mm (red)	Bag/100 pcs
8109432	8 mm (blue)	Bag/100 pcs
8109433	9 mm (green)	Bag/100 pcs
8109434	10 mm (yellow)	Bag/100 pcs
8109435	11 mm (red)	Bag/100 pcs
8109436	12 mm (blue)	Bag/100 pcs
8109437	13 mm (green)	Bag/100 pcs
8109438	14 mm (yellow)	Bag/100 pcs
8109439	16 mm (red)	Bag/100 pcs
8109440	18.5 mm (blue)	Bag/100 pcs



Impedance ear tips

Sanibel™ ATE umbrella ear tips

For use with old MT10™ and the following systems with BLACK probe tip: AT235™, AA220™, MI22™, MI23™, MI24™, MI26™, MI34™, MI44™ and old SPL60 (4 mm).

Part no.	Description	Unit
8107617	10 mm (green)	Bag/100 pcs
8107618	13 mm (green)	Bag/100 pcs
8107619	16 mm (green)	Bag/100 pcs
8107620	19 mm (blue)	Bag/100 pcs
8107621	22 mm (blue)	Bag/100 pcs



Impedance ear tips

Sanibel™ AZE mushroom ear tips

For use with old MT10™ and the following systems with BLACK probe tip: AT235™, AA220™, MI22™, MI23™, MI24™, MI26™, MI34™, MI44™ and old SPL60 (4 mm).

Part no.	Description	Unit
8106189	6 mm (violet)	Bag/25 pcs
8107622	6 mm (violet)	Bag/100 pcs
8106190	7 mm (yellow)	Bag/25 pcs
8107623	7 mm (yellow)	Bag/100 pcs
8106191	9 mm (green)	Bag/25 pcs
8107624	9 mm (green)	Bag/100 pcs
8106192	10 mm (blue)	Bag/25 pcs
8107625	10 mm (blue)	Bag/100 pcs
8106193	11 mm (white)	Bag/25 pcs
8107626	11 mm (white)	Bag/100 pcs
8106194	13 mm (yellow)	Bag/25 pcs
8107627	13 mm (yellow)	Bag/100 pcs
8106195	15 mm (green)	Bag/25 pcs
8107628	15 mm (green)	Bag/100 pcs
8106196	18 mm (blue)	Bag/25 pcs
8107629	18 mm (blue)	Bag/100 pcs



Ear tips for OAE & ABR screening

Sanibel™ ANA ear tips

For use with Echo-Screen®, AccuScreen® Classic and Eclipse™ with pink probe tip.

Part no.	Description	Unit
8106363	4.0 mm (blue)	Bag/50 pcs
8106362	4.5 mm (red)	Bag/50 pcs
8106364	5.0 mm (green)	Bag/50 pcs
8106366	Infant tree (yellow)	Bag/50 pcs
8106365	Mushroom (green)	Bag/50 pcs



Ear tips for OAE & ABR screening

Sanibel™ AOT ear tips

For use with AccuScreen® (touch) and Sentiero™ with small probe tip.

Part no.	Description	Unit
8106367	3.7 mm (yellow)	Bag/50 pcs
8106368	4.0 mm (light green)	Bag/50 pcs
8106369	4.5 mm (light pink)	Bag/50 pcs
8106370	5.0 mm (light blue)	Bag/50 pcs
8106371	4-7 mm (transparent)	Bag/50 pcs



Ear tips for OAE screeners

For use with AuDX®, ABaer®, Navigator® Pro, and Echo-Screen® III.

Part no.	Description	Unit
8107885	AuDX Preemie 3 mm	Bag/100 pcs
8515505	AuDX Preemie 3 mm	Bag/20 pcs
8107889	AuDX Tree 3.5-5mm	Bag/100 pcs
8515506	AuDX Tree 3.5-5mm	Bag/20 pcs
8107887	AuDX Infant 6 mm	Bag/100 pcs
8515507	AuDX Infant 6 mm	Bag/20 pcs
8102681	Adult OAE foam	Bag/100 pcs
8515436	Adult OAE foam	Bag/20 pcs
8102682	Pediatric OAE foam	Bag/100 pcs
8515435	Pediatric OAE foam	Bag/20 pcs
8103123	Jumbo OAE foam	Bag/100 pcs
8515438	Jumbo OAE foam	Bag/20 pcs



Grason & Associates OT select ear tips w/two probe holes

For use with Otodynamics® ECP™, UGS™, SNS™ & SGS™ probes.

Part no.	Description	Unit
8102632	3.5 mm (blue)	Bag/100 pcs
8102633	4.5 mm (blue)	Bag/100 pcs
8102634	5.5 mm (blue)	Bag/100 pcs
8102635	6.5 mm (blue)	Bag/100 pcs
8102636	4.8 mm (blue)	Bag/100 pcs
8102637	6 mm (green)	Bag/100 pcs
8102638	8 mm (red)	Bag/100 pcs
8102639	10 mm (blue)	Bag/100 pcs
8102640	12 mm (green)	Bag/100 pcs
8102641	14 mm (red)	Bag/100 pcs
8102642	16 mm (blue)	Bag/100 pcs



Grason & Associates OT select ear tips w/triangular hole

For use with Otodynamics® UGD™, SND™ & SGD™ probes.

Part no.	Description	Unit
8102644	4.8 mm (blue)	Bag/100 pcs
8102645	6 mm (green)	Bag/100 pcs
8102646	8 mm (red)	Bag/100 pcs
8102647	10 mm (blue)	Bag/100 pcs
8102648	12 mm (green)	Bag/100 pcs
8102649	14 mm (red)	Bag/100 pcs
8102650	16 mm (blue)	Bag/100 pcs



Grason & Associates MO ear tips

For MADSEN® OTOflex™.

Part no.	Description	Unit
8102653	3-6 mm (yellow)	Bag/100 pcs
8102654	6 mm (red)	Bag/100 pcs
8102655	7 mm (blue)	Bag/100 pcs
8102656	8 mm (green)	Bag/100 pcs
8102657	9 mm (yellow)	Bag/100 pcs
8102658	10 mm (red)	Bag/100 pcs
8102659	11 mm (blue)	Bag/100 pcs
8102660	12 mm (green)	Bag/100 pcs
8102661	13 mm (yellow)	Bag/100 pcs
8102662	14 mm (red)	Bag/100 pcs
8102663	16 mm (blue)	Bag/100 pcs
8102664	18 mm (green)	Bag/100 pcs



Save time pr. patient, earn more money and help more patients!

Use a fresh ear tip pr. patient and be sure to reach a seal right away, when inserting the ear tip. A fresh ear tip contains the optimal conditions for reaching a seal when inserting it into the patients' ear, whereas a multi-used one often will cause you trouble and delay you in

your process. By allowing yourself to plan less time pr. patient for testing, you can see more patients pr. day. That is much more worth than an ear tip, both for your clinic and your patients.



We know that, first and foremost, great quality is at the top of your list, and you'll find it in each and every purchase you make from Sanibel Supply.

Bringing you the best is truly what drives us day after day.

If we can ever be of assistance, don't hesitate to let us know, contact@sanibelsupply.com


Jakob Skovgaard
 Business Unit Manager

Probe tips



IMP

Probe tip for Sanibel™ ADI™ ear tips

Fits Otowave and new MT10™.

Part no.	Unit
8002592	1 pc



ABR IMP OAE

Probe tip for Sanibel™ ADI™ ear tips

Fits Titan™ and easyTymp™.

Part no.	Unit
8001701	1 pc



IMP

Probe tip kit

Fits AT235™ or newer, AA222™, Sera™, Tympstar Pro™ Novus™, easyScreen™, and touchTymp™. Includes probe tip and gasket.

Part no.	Unit
8107762	Set/2 pcs



ABR OAE

Sanibel™ probe tubes for ADI™ ear tips

Fits Corti™, ERO•SCAN™ and the new OtoRead™.

Part no.	Unit
8104159	Box/100 pcs



ABR OAE

Probe Tube Removal Tool

Fits Corti™, ERO•SCAN™ and the new OtoRead™.

Part no.	Unit
8513131	Bag/1 pc



ABR OAE

Probe Tube Removal Tool kit

Fits Corti™, ERO•SCAN™ and the new OtoRead™.

Part no.	Unit
8513132	Bag/100 pcs + 1 pc



OAE

Probe tip kit for OAE Ø3.3 mm ear tips

Fits all Eclipse™ – use with ADI™ ear tips. Requires no re-calibration.

Part no.	Unit
8012960	Bag/8 pcs



IMP

Sanibel™ probe tip adaptor kit

Converts AT235™, AT235h™, AA220™, AA222™(blue/new), MI24™, MI26™, MI34, MI44™ for use with ADI™ ear tips. Includes probe tip, washer, damper. No recalibration required.

Part no.	Unit
8500430	Set/3 pcs



OAE

Remote probe tip replacements for ERO•SCAN™

Fits ERO•Scan™.

Part no.	Unit
8102340	Bag/4 pcs





Internal probe tip replacements for ERO•SCAN™

Fits ERO•SCAN™.

Part no.	Unit
8102339	Bag/4 pcs



ERI0D probe tips for ERO•SCAN™ Pro and OtoRead™

It is recommended to clean the probe tip after every patient and to change the probe tip after 20 tests (40 ears). No calibration necessary.

Part no.	Unit
8012836	Bag/4 pcs



Internal probe tip replacements for ERO•SCAN™ Pro

Fits ERO•Scan™ Pro.

Part no.	Unit
8110364	Bag/4 pcs



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags

3A E-A-RLINK™ foam ear tips

Complies to all insert phones.


Part no.	Description	Unit
8500090	Adult/Medium size 13.7 mm	Bag/50 pcs
8106226	Adult/Medium size 13.7 mm	Bag/100 pcs



E-A-RLINK™ infant ear tips

For use with infants or in very small ear canals.

Part no.	Description	Unit
8500300	3.5 mm	Bag/20 pcs
8500310	4.0 mm	Bag/20 pcs



Sanibel™ IP30™ tube connectors

For use with IP30™ and E-A-Rtone® 3A.

Part no.	Unit
8107391	Bag/10 pcs

Sanibel™ Probe Tube Removal Tool


- A = Probe
- B = Probe Tube Installer
- C = Attachment to probe cable
- D = Tube Remover

To Remove a Probe Tube

- Place the front hole of the Probe Tube Tool over the end of the probe, as shown. The probe should be seated against the face of the tool and snap in place.
- Squeeze the tool closed and hold it.
- Twist the tool a couple of times while holding the tool closed and pulling away from the probe.
- The probe tube will pull out from the probe. Discard the probe tube.

To Install a New Probe Tube

- Place a new tube into the hole on the top of the tool.
- Align the probe with the tube.
- Push the probe onto the probe tube until the probe contacts the face of the probe tube tool and a snapping sound is heard.
- Slide the tool off leaving the tube in place in the probe.



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags

3B E-A-RLINK™ foam ear tips

Complies to all insert phones.


Part no.	Description	Unit
8500110	Pediatric/Small size 9.7 mm	Bag/50 pcs
8106227	Pediatric/Small size 9.7 mm	Bag/100 pcs



Sanibel™ ADI ear tip adaptor

Makes Sanibel™ ADI ear tips compatible with 3A insert phones.

Part no.	Unit
8500580	Bag/6 pcs



Tube nipples for E-A-Rtone® 3A

For use with sound tubes 3A.

Part no.	Unit
8500140	Bag/10 pcs



3C E-A-RLINK™ foam ear tips

Complies to all insert phones.

Part no.	Description	Unit
8500100	Jumbo/Large size 17.8 mm	Bag/24 pcs



Sound tubes for E-A-Rtone® 3A ABR 50 ohm

Includes adaptor nipples, 1 pc blue and 1 pc red 50 ohm.

Part no.	Unit
8500720	Set/2 pcs



Sanibel™ sound tubes and Sanibel™ ADI™ ear tip adaptor kit

Fits E-A-Rtone® 3A inserts and IP30™ insert phones. Includes ADI™ ear tip adaptors, 1 pc blue and 1 pc red 50 ohm.

Part no.	Unit
8100590	Set/2 pcs



ABR
Sound tubes for E-A-Rtone® 3A 10 ohm
 Includes adaptor nipples.

Part no.	Unit
8107396	Bag/4 pcs



ABR
Sound tubes for E-A-Rtone® 3A 10 ohm
 Includes adaptor nipples.

Part no.	Unit
8103605	Bag/4 pcs



ABR
Sanibel™ IP30 and E-A-Rtone® sound tubes and nipples
 Fit all IP30™ and E-A-Rtone® versions. Includes tube connectors, 1 pc blue and 1 pc red 50 ohm.

Part no.	Unit
8501072	Set/2 pcs



ABR
E-A-Rtone® replacement fastener clip
 Clip for fastening insert phone on clothes.

Part no.	Unit
8000744	1 pc



ABR
Sanibel™ fastener kit for IP30™
 Includes velcro neckstrap and clip.

Part no.	Unit
8507410	Set



ABR OAE
Sanibel™ ear tip adaptor cleaning tool
 Small brush for cleaning the Sanibel™ probe tip adaptors before disinfection.

Part no.	Unit
8500830	Bag/10 pcs

100 tests



ABR
ABR consumable kit
 Includes electrodes, insert phone tips, skin scrub and gaze.

Part no.	Description	Unit
8109670	Adult 3A	Kit/100 tests
8109671	Pediatric 3B	Kit/100 tests

2

ABR
Sanibel™ Infant EarCup™ with snap or tab electrodes

Reduce test time! Test on both ears in one screening. Compatible with Titan™, Eclipse™, AccuScreen® (1st generation), MB II Classic and Echo-Screen®. Use only with equipment calibrated for over-ear style disposables. Disposable over-the-ear earphones are high quality, low cost, and ideal for bilateral newborn hearing screenings. 1 set consists of 2 Infant EarCups™, 3 snap or tab electrodes and 1 alcohol pad.

Part no.	Description	Unit
8049600	Tab	Box/20 sets
8049700	Snap	Box/20 sets



Exchange the electrodes with your preferred version!
CONTACT US FOR SUPPORT



ABR
Sanibel™ Infant EarCup™ adaptor set
 Connects insert phones with Infant EarCups™. NB: Do not contain filters.

Part no.	Unit
8500450	Set/2 pcs



ABR
Sanibel™ sound tubes and Sanibel™ Infant EarCup™ adaptor kit
 Includes adaptors, 1 blue and 1 red tube.

Part no.	Unit
8100577	Set/2 pcs





Use Sanibel™ probe tubes for any equipment

We often hear from the industry that only the manufacturers' specific probe tubes should be used with their specific equipment. However, a recent comparative study by Sanibel Supply® disproves this myth.

In this study, our in-house audiologists and engineers compared the effects of different probe tubes on an Affinity^{2.0} system. They compared AudioScan tubes (red), Otometrics tubes (blue), Sanibel™ tubes (orange) and Sanibel™ tubes (grey).

This involved running each tube via the 'tube calibration' process and then comparing the following REUG curve when the tube had been left in the calibration position.

The investigation showed no significant difference in the REUG curve and all probe tubes achieved the same result when the calibration process had been run correctly.

The "REM effect"
– better fitting & fewer follow-ups!

- Be sure of an individual and ideal fitting
- Assurance of your fitting gives you confidence in your work
- Confidence increases your credibility and client satisfaction
- Client satisfaction means fewer follow-up appointments



HIF
Sanibel™ REM probe tubes, type A
Universal fit, length 76 mm, beveled 45° angle tip and black marker ring.

Part no.	Unit
8104116	Tray/50 pcs



HIF
Sanibel™ REM probe tubes, type B
Cost efficient, length 80 mm and marker ring.

Part no.	Unit
8500360	Bag/40 pcs



HIF
Silicone tubes for SPL60™
For use with Affinity™ REM probe.

Part no.	Unit
8500700	Bag/50 pcs



Large soft electrode size! 60 pcs

ABR
Sanibel™ snap electrodes
Pre-gelled with a very gentle hydrogel, perfect for ABR testing on infants. Size 2.2 x 3 cm. Item no. 8107137 is packed in 20 smaller pouches of 3 electrodes.

Part no.	Unit
8107137	Box/60 pcs (20x3 pcs)
8109096	Bag/60 pcs (10x6 pcs)



Large soft electrode size! 60 pcs

ABR
Sanibel™ tab electrodes
Pre-gelled cloth electrode with a very gentle hydrogel, perfect for ABR testing on infants. Size 2.2 x 3 cm. Packed in 20 smaller pouches of 3 electrodes.

Part no.	Unit
8107410	Box/60 pcs (20x3 pcs)



ABR
Sanibel™ wetgel snap electrodes
Pre-gelled round foam ECG electrode with button connector. Packaged in convenient 12x5 pcs re-sealable pouches.

Part no.	Unit
8108563	Bag/60 pcs (12x5 pcs)



ABR
Sanibel™ vinyl tab electrode
For alligator clip cable. Solid gel. Size 2.5 x 2.2 cm. Packaged in convenient re-sealable pouches.

Part no.	Unit
8107426	Bag/100 pcs (10x10 pcs)



ABR
Sanibel™ TM electrode
For ECochG. Use with cable 8500650.

Part no.	Unit
8011489	Box/10 pcs



ABR
Sanibel™ prewired electrodes
Prewired electrodes with a small footprint. Very good for testing infants. Solid gel. Size 2.2 x 2.2 cm. 100 pouches of 3 pcs.

Part no.	Unit
8107425	Box/300 pcs



ABR
E-A-RLINK™ gold foil tip trodes
For ECochG. Use with cable 8011343.

Part no.	Description	Unit
8500370	Adult 13 mm	Bag/20 pcs
8500440	Pediatric 10 mm	Bag/20 pcs



ABR
Vermed Economy line 3-pak Solid gel snap Electrode
Solid Gel type snap electrode suitable for newborn screening and Diagnostic follow-up. Latex free. Intended use areas: Newborn, Diagnostic follow-up all ages.

Part no.	Unit
8516252	Box/1200 pcs (20x60 pcs)

Remember to prep your patients for great results



60 pcs

ABR

Vermed Wetgel Snap Electrode

Pre-gelled hydrogel electrode with Ag/AgCl sensor. Very suitable for newborn screening or diagnostic follow up. Packaged 3x20 pcs. per pouch.
Intended use areas: Newborn, Diagnostic follow-up, ABR, C-Vemb.

Part no.	Unit
8516221	Bag/60 pcs (3x20 pcs)



1000 pcs

ABR

Vermed Tab Electrode

Rectangular tab electrode for newborn screening. 12,7x22,9 mm. Conductive Adhesive Hydrogel. 10x10 pcs. per pouch.
Intended use areas: Newborn, O-Vemb.

Part no.	Unit
8516225	Bag/1000 pcs (100 pcs per pouch)



60 pcs

ABR

Vermed Hydrogel Snap Electrode

Large soft electrode size making it very suitable for newborn screening and ABR measurements. Cloth material makes it convenient to wear. Size: 22.2.
Intended use areas: Newborn, ABR.

Part no.	Unit
8107137	Bag/60 pcs (3x20 pcs)



100 pcs

ABR

Ambu® 0415M tab electrode

Medium size tab electrode with solid gel for 12-lead ECG. Presented conveniently with 10 pieces to a card. Fish tail tab design reduces electrode lift-offs. Suitable for both adults and children.
Intended use areas: ABR Diagnostics of all kinds with Tab electrodes.

Part no.	Unit
8516138	Bag/100 pcs (10x10 pcs)



1000 pcs

ABR

Ambu® 2300 Blue Sensor Tab electrode

The medium size Ambu BlueSensor 2300 electrode features a skin-friendly solid gel that provides optimal signal quality recordings during 12-lead diagnostic resting ECG application. Very suitable for ABR diagnostics of all kinds.
Intended use areas: ABR Diagnostics of all kinds.

Part no.	Unit
8516139	Box/1000 pcs (10x100 pcs)



25 pcs

ABR

Ambu® Blue Sensor™ N-00-S

Foam ECG electrode for use with button cables. Wetgel.

Part no.	Unit
8013029	Bag/25 pcs



60 pcs

ABR

Vermed Hydrogel Snap Electrode

Smaller than the rectangular version, but otherwise similar. Good soft electrode size making it very suitable for newborn screening and ABR measurements. Cloth material makes it convenient to wear.
Intended use areas: Newborn, ABR.

Part no.	Unit
8102022	Bag/60 pcs (3x20 pcs)



1000 pcs

ABR

Vermed Tab Electrode

Rectangular tab electrode for newborn screening. 22.9x34 mm. Conductive Adhesive Hydrogel. 10x10 pcs. per pouch.
Intended use areas: Newborn, O-Vemb.

Part no.	Unit
8516227	Bag/1000 pcs (100 pcs per pouch)



60 pcs

ABR

Vermed Wet Gel foam

Electrode with an Ag/AgCl sensor. 10x6 pcs per pouch. Diameter 38.1 mm

Part no.	Unit
8108563	Bag/60 pcs



25 pcs

ABR

Ambu® Neuroline 720

Ambu Neuroline 720 electrodes are pre-gelled and self adhesive. They are ideal for use in sleep registration where the strong, yet skin-friendly, adhesion ensures that the electrode stays in place for the entire examination, even if the patient sweats.
Intended use areas: ABR Diagnostics of all kinds.

Part no.	Unit
8500060	Bag/25 pcs



50 pcs

ABR

Ambu® Blue Sensor™ P-00-S

High quality Ag/AgCl sensor – ensures excellent trace quality.

Part no.	Unit
8504783	Bag/50 pcs



ABR

Ambu® Neuroline 720 - with 100 cm cable and k-connector

Ambu Neuroline 720 electrodes are pre-gelled and self adhesive. They are ideal for use in sleep registration where the strong, yet skin-friendly, adhesion ensures that the electrode stays in place for the entire examination, even if the patient sweats.
Intended use areas: ABR Diagnostics of all kinds.

Part no.	Unit
8503244	Bag/10 pcs



300 pcs

ABR

Vermed Prewired neonatal ecg

electrode for use in pediatric and neonatal settings. 3 electrodes per pouch. square shape 22.2x22.2 mm

Part no.	Unit
8107425	Bag/300 pcs (3x100 set)




1000 pcs

ABR

Vermed Economy line Snap Wet Gel Electrode

Non-woven cloth electrode with conductive wet gel. Passed ISO 10993 and is Latex free.
Intended use areas: Newborn, Diagnostic follow-up all ages.

Part no.	Unit
8516251	Box/1000 pcs (20x50 pcs)



60 pcs

Vermed Allround ECG electrode

Provides very good impedance and it well suited for both newborns and adults. Can be used both as snap and tab electrode with alligator cables. Pouches come in 3 pcs, value packs available in 600 pcs. per order. Please contact for order.

Part no.	Unit
8516906	Bag/60 pcs (5x12 pcs)



25 pcs

ABR

Ambu® Blue Sensor™ N-10-A

Foam ECG electrode with 10 cm wire and 4 mm connector, for use with universal electrode cable. Wetgel.

Part no.	Unit
8500670	Bag/25 pcs



New selection of surface electrodes for electrophysiological measurements.

Our selection of electrodes has been widened to ensure your needs are met at an attractive price.

The electrodes are compatible with all electrophysiological brands.

However, Sanibel Supply surface electrodes are not suitable for measuring the heart rate, they are specific tailored for auditory brainstem responses.

Measuring auditory brainstem response requires very sensitive sensors, whereas less sensitive sensors can measure the heart rate.

As such, various surface electrodes are suited for various measurements.

Conflicting surface electrodes and measurements will lead to misleading results, and ultimately a false diagnosis.

Surface electrodes tailored to your need

At Sanibel Supply® we believe that it is important to choose the right surface electrode for each specific recording and patient type.

Wetgel or hydrogel electrodes?

A wetgel electrode consists of water, thickener, ionic salts, surfactants and bactericides or fungicides. Wetgel migrates easily and rapidly onto the skin surface, reducing impedance levels, which in turn will generally result in lower noise in the recording.

A hydrogel surface electrode is less adhesive when compared to wetgel, making it particularly suitable for newborns.

Chloride salt concentration

The higher the chloride salt content, the more conductive the electrode. Higher salt content (10% and above), pre-gelled surface electrodes are useful for making fast, high quality measurements of biopotentials. High conductivity electrodes generally have reduced artifact, due to the lower impedance generated between electrode and skin surface. As the chloride salt content drops, there is a corresponding reduction in electrode conductivity.

If the chloride content is below 10%, the electrode can then be employed for longer recording times, with a reduced risk of skin irritation. As a rule of thumb, wetgel has a chloride content of 10% and hydrogel 4%.

It is important to use electrodes with similar chloride salt content when testing.

The impedance of the electrode/skin junction is highly dependent on the electrolyte type and the chloride salt concentration. For example, a hydrogel electrode with 4% chloride concentration will have approximately 10x higher impedance than a wet liquid gel electrode with 10% chloride concentration, after first application to the skin.

Adhesion

An advantage of a conductive adhesive is that the active surface can be as large as the electrode itself because no additional adhesive strip is needed. The Sanibel disposable electrodes come with a range of adhesive qualities. All electrodes are designed to adhere well to skin surfaces, but


the adhesion characteristics can be grouped according to the application:

- Strong adhesive
- Moderate adhesive

Strong adhesion is not recommended for newborns.

Skin preparation

For optimum electrode to skin conductivity, the skin should be lightly brushed with a gentle abrasive gel/wipe. **An alcohol wipe is not recommended to improve conductivity as this will only dry out the skin surface.**




Choose the right electrode for optimal performance

Sanibel Supply® offers various kinds of surface electrodes for ECG as well as EEG and ABR. Surface electrodes vary in size, adhesion, connectivity and frequency (test category).

We test our surface electrodes in a controlled setup to ensure they are approved for the test category in question.

The frequency plot to the left shows the electrical stimuli sent directly through the system and compared to the surface electrode (here ABR). The two lines are similar and demonstrate that the surface electrode can be used for ABR testing.

At Sanibel Supply®, we always test and validate our surface electrodes to ensure optimal performance for your testing!

Item no.	Surface electrode	Adhesion electrode backing / Dimension (part on skin) / Material hardness / Electrode cable	Chloride conc. - Gel type	Category / Patient age	Special conditions
8107410	 Sanibel™ tab electrodes Packed in 20 smaller pouches of 3 electrodes (total of 60 electrodes).	Moderate / 3 x 2.2 cm / Soft tissue / Tab cables (ETSE tab surface cable kit #8011352)	Solid gel / Hydrogel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Low chloride content, less adhesive. Especially suitable for newborns.
8107137	 Sanibel™ snap electrodes Pouch of 3 electrodes bundled in a box with a total of 60 electrodes.	Moderate / 3 x 2.1 cm / Soft tissue / Snap cable (ETB surface electrode kit #8011346)	Solid gel / Hydrogel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Low chloride content, less adhesive. Especially suitable for newborns.
8107425	 Sanibel™ prewired electrodes 3 electrodes in each pouch with a total of 300 electrodes.	Moderate / 2 x 2 cm / Soft tissue / Prewired	Solid gel / Hydrogel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Low chloride content, less adhesive. Especially suitable for newborns. Electrode cables are disposable, reducing cross contamination. Small soft electrode.
8108563	 Sanibel™ wetgel snap electrodes 12x5 pcs re-sealable pouches (total of 60 electrodes).	Strong / Ø 1 inch, 2.55 cm / Soft plastic / Snap cable (ETB surface electrode kit #8011346)	3% - Wetgel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Strong adhesion. Low chloride content (not recommended for newborns under 1 month of age).
8109096	 Sanibel™ snap electrodes 60 electrodes in one resealable pouch. Repacked in 10 smaller bags of 6 pcs.	Moderate / 3.5 x 2.3 cm / Soft tissue / Snap cable (ETB surface electrode kit #8011346)	Solid gel / Hydrogel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Low chloride content, less adhesive. Especially suitable for newborns. Large soft electrode.
8107426	 Sanibel™ vinyl tab electrodes 100 electrodes in one resealable pouch.	Moderate / 2.5 x 2.2 cm / Hard plastic / Tab cables (ETSE tab surface cable kit # 8011352)	Solid gel / Hydrogel		Low cost tab electrode.
8500060	 Ambu® Neuroline™ 720 side snap 25 electrodes in one pouch.	Strong / 2.8 x 2 cm / Soft plastic / Snap cable (ETB surface electrode kit #8011346)	Wetgel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Strong adhesion. High chloride content (not recommended for newborns under 1 month of age).





Sanibel™ ECoHG starter kit
For Eclipse™ and Bio-logic® AEP system. Includes cable, gel and 2 pcs TM electrodes.

Part no.	Unit
8013209	Set/4 pcs



Sanibel™ ECoHG cable
For Eclipse™ and Bio-logic® AEP system. 1.5 mm/.060" female connector. For use with Sanibel™ TM electrodes.

Part no.	Unit
8500650	1 pc



Sanibel™ EP electrode jumper cable
Allows two of the preamplifier inputs to be connected, e.g. for linking mastoid electrodes for cortical testing.

Part no.	Unit
8107872	1 pc



Tip trode cables
For use with TTE25 gold foil tip trodes. 1.5 mm with touch proof connectors and plated clamp clips. 50 cm long.

Part no.	Unit
8011343	Set/2 pcs



NuPrep™ skin prep gel, 114 g
Efficient, yet gentle skin abrasive. Chloride and acetone free. Pack of 1x4 oz/114 g tube or 3x4 oz/114 g tube.

Part no.	Unit
8500010	1 tube
8500320	Box/3 tubes




Lemon Prep™ abrasive lotion, 114 g
Pack of 3x4 oz/114 g tubes.

Part no.	Unit
8500770	Box/3 tubes

After use ...


After use, cables should be wiped with a soft cloth that has been dipped in a mild detergent solution or alcohol.

Contact between the cleaning solution and the connectors should be avoided.



Disc electrodes
For use with Ten20™ electrode paste. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable.

Part no.	Unit
8011498	Set/3 pcs RWB
8011345	Set/5 pcs RWBYK



Universal electrode cable
Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable. Non-tangling.

Part no.	Unit
8011353	Set/3 pcs RWB
8011344	Set/5 pcs RWBYK



Lectron® II electrode gel, 250 ml
Highly conductive, chloride-free electrode gel. Recommended for use with Sanibel™ TM electrodes for ECoHG.

Part no.	Unit
8102626	Tube/250 ml




SignaGel®, 250 g
Hypoallergenic, bacteriostatic conductive saline electrode gel.

Part no.	Unit
8102628	Tube/250 g



Gold-plated cup electrodes
Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). Ø 10 mm gold-plated. 50 cm cable.

Part no.	Unit
8107418	Set/3 pcs RWB
8107145	Set/5 pcs RWBYK



Alligator clip cable
For use with tab electrodes. 1.5mm/.060" female connector. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K).

Part no.	Unit
8011412	Set/3 pcs RWB
8500900	Set/3 pcs WYK
8011352	Set/5 pcs RWBYK



Pinch style cable
For use with snap electrodes. 1.5mm/.060" female connector. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K).

Part no.	Unit
8500390	Set/3 pcs RWB
8500400	Set/3 pcs WYK
8500410	Set/5 pcs RWBYK



Standard button cable
For use with snap electrodes. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable. Non-tangling.

Part no.	Unit
8011354	Set/3 pcs RWB
8011346	Set/5 pcs RWBYK



Spectra 360® electrode gel, 250 g
Salt- & chloride-free conductive electrode gel suitable for long-term applications.

Part no.	Unit
8500740	Tube/250 g



Ten20™ conductive EEG paste, 114 g/228 g
Adhesive gel for disc electrodes. Available in pack of 3x4 oz/114 g tubes, 3x8 oz/228 g jar or 1x8 oz/228 g jar.

Part no.	Unit
8012924	Box/3 tubes
8500070	Box/1 jar
8500330	Box/3 jars



ABR
Everi™ abrasive paste, 160 g
 Reduces the skin impedance during EEG, EP, ECG and Holter Stress test examination.

Part no.	Unit
8012952	Bottle/160 g



ABR
Sonavelle® electrode gel, 250 ml
 Formaldehyde-free and glycerol-free.

Part no.	Unit
8500380	Bottle/250 ml



ABR
3M™ Red Dot™ trace prep, 5 m
 One-step skin preparation. 18 mm x 5 m / 0.70" x 196".

Part no.	Unit
8500030	Dispenser/5 m



12 rolls

ABR
3M™ Transpore™ Surgical tape, 9.14 m
 Transparent, latex-free, water-resistant and porous for breathability. 25 mm x 9.14 m / 1" x 10 LY.

Part no.	Unit
8500750	Box/12 rolls



12 rolls

ABR
3M™ Micropore™ Surgical tape, 9.14 m
 Latex-free and hypoallergenic. Holds well on damp skin. 25 mm x 9.14 m / 1" x 10 LY.

Part no.	Unit
8500760	Box/12 rolls



10 packs

ABR OAE
Q-tips
 Cotton tipped applicators for patient preparation. Non-sterile. Each pack contains 100 pcs.

Part no.	Unit
8500810	Box/10 packs



ABR OAE
Gauze swabs
 Non-sterile, 7.5 x 7.5 cm and 8 layered.

Part no.	Unit
8500040	Pack/100 pcs.



ABR OAE
B Braun® alcohol pads
 70% isopropyl alcohol.

Part no.	Unit
8500050	Box/100 pcs.



500 pcs



ABR AUD
3" Disposable ear cushion cover
 Replaces disinfectants. Headphone Covers for classic audiometric headphones.

Part no.	Unit
8107419	Box/500 pcs



100 pcs



ABR AUD
4.5" Disposable ear cushion cover
 Replaces disinfectants. Headphone Covers for noise enclosure audiometric headphones. For HDA200/DD450/DD65v2/ Peltor headsets.

Part no.	Unit
8513404	Bag/100 pcs



Keyboard protective foil
 Place a protective foil over your keyboard to minimize contamination with grime, bacteria and germs, or liquid.

Part no.	Unit
8500460	1 pc



ABR IMP OAE
Probe cleaning tool
 For Interacoustics Titan™/AT235™/AA222™/ Sera™/MT10™/Amplivox Otowaves/MAICO Mi XX + touchTymp™ and GSI TympStar™ Pro.

Part no.	Unit
8011351	1 pc



Disinfectant detergent, 1000 ml
 Alcohol-free, ready to use solution for a quick and extensive disinfection of non-invasive medical devices.

Part no.	Unit
8501706	Bottle, 1000 ml



ABR IMP OAE
Bridge & Implant floss
 Used to clean probe tips.

Part no.	Unit
8003566	Box/30 pcs



How to clean probe tips!



AUD
USB cables A to B

Part no.	Description	Unit
8011241	2 m / 6.56 feet	1 pc
8011243	5 m / 16.40 feet	1 pc



AUD
TDH39™/DD45™ replacement cable, symmetric
 For TDH39™ and DD45™. 30° jack plug.

Part no.	Description	Unit
8010822	2 m	1 pc
8507204	5 m	1 pc




Medical grade silicone rubber
Biocompatible
AUD
Sanibel™ TDH39/DD45 earphone cushion red & blue, model 51
 For TDH39, 49, 50 & DD45 earphones.

Part no.	Unit
8107611	1 Set





AUD
RadioEar Foldable audiometric headband
 Black, foldable.

Part no.	Unit
8010882	1 set



AUD
RadioEar Headband with 2 yokes
 Black.

Part no.	Unit
8010840	1 set



AUD
Peltor® headset/earmuff replacement hygiene kit
 For DD450, HDA200 and DD65.

Part no.	Unit
8500270	1 set



AUD **IMP**
TDH39™/DD45™ contra headset replacement cable
 For TDH39™ and DD45™. Single mono. 30° jack plug.

Part no.	Unit
8010826	1 pc



AUD
TDH39™/DD45™ replacement cable, non-symmetric
 For TDH39™ and DD45™. Non-symmetric mono. 30° jack plug.

Part no.	Unit
8010823	1 pc




Medical grade silicone rubber
Biocompatible
AUD
DD 51 cushion
 For DD45™ and TDH39™ headsets.

Part no.	Unit
8010894	1 pc





AUD
RadioEar DD45 Audiometric headset
 Headband: HB7. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8010954	1 set



AUD
RadioEar DD45 Adult audiometric headset
 Headband: HBA, foldable. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8106340	1 set



AUD
HDA300 replacement cable, symmetric
 For HDA300 and DD65 headphones. 2 m cable and 6.3 mm/ 30° jack plug.

Part no.	Unit
8004172	1 pc



AUD
B71 Bone conductor replacement cable
 2 m cable and 6.3 mm/ 30° jack plug.

Part no.	Unit
8011094	1 pc



AUD
IP30 Contra Insert earphone cable
 Cable length: 2 m. Plug type: 30°, 1x6.3 mm gold-plated mono jack.

Part no.	Unit
8107358	1 pc



AUD
RadioEar DD45 Adult audiometric headset
 Headband: HB3045. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8105268	1 set





AUD
RadioEar DD65 v2 Passive attenuation audiometric headset

Cable length: 2 m. Plug type: Straight, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8510117	1 set



AUD
RadioEar B71W Bone transducer headset

Ohm: 10 ohm. Cable length: 79"/200 cm molded cord. Plug type: 30°, 1x6.3 mm gold-plated mono jack. Headband: P3333

Part no.	Unit
8510370	1 pc



AUD
RadioEar IP30 Insert earphone cable

Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8004432	1 pc



AUD
Foam pad for contra headset

Compatible with e.g. TDH39c™ and DD45c™.

Part no.	Unit
8507292	Bag/10 pcs



AUD HIF ABR
RadioEar replacement pads (3288-1)

For RadioEar Spring Steel headband.

Part no.	Unit
8504777	Bag/10 pcs



AUD
Ear cushion HZP 09 pair

Ear cushion foam (1 pair) suitable for: PC 20, PC 25, PC 30, PC 31, PC 35, PC 36, PC 130, PC 131, PC 135, PC 136, PC 140, PC 141, PC 145, PC 146.

Part no.	Unit
8104416	1 pair



AUD
Headband cushion

For RadioEar HBA foldable audiometric headband.

Part no.	Unit
8010886	1 pc



AUD
Headband cushion

For Sennheiser HDA200, RadioEar DD450, DD65v2 and Peltor audiometric headbands.

Part no.	Unit
8000770	1 pc



AUD HIF ABR
RadioEar Spring Steel headband (P3333)

Part no.	Unit
8507213	1 pc



AUD
RadioEar DD45 transducer

Audiometric earphones engineered to match the performance and function of the most popular earphones available for audiological testing. 10 ohm receiver. Easy installation. Reliable. RoHS compatible.

Part no.	Unit
8507268	1 set



AUD
RadioEar DD45 Contra headband

Tension: 4,5 N ±0.5N

Part no.	Unit
8507211	1 set



AUD HIF ABR
RadioEar AMBAND bone conductor headband (P1188)

For adult, children and infants.

Part no.	Unit
8507214	Bag/10 pcs





24 pcs

Disposable VNG goggle foam pads

For Interacoustics® VNG mask.

Part no.	Unit
8011527	Box/24 pcs



34 pcs

Disposable ear specula

For AirFx™.

Part no.	Description	Unit
8013013	2.75 mm	Bag/34 pcs
8013014	4.25 mm	Bag/34 pcs



50 pcs

Water irrigator tips

For use with Interacoustics® Aqua Stim™.

Part no.	Unit
8013018	Bag/50 pcs



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 rolls

OAE

Thermal paper for ERO-SCAN™/OtoRead™, 25 m

Width: 57 mm
Inner diameter: 12 mm
Outer diameter: 45 mm

Part no.	Unit
8500200	1 roll/25 m



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 rolls

IMP

TPR10™ thermal paper for MTP10™, 25 m

Width: 79 mm
Inner diameter: 12 mm
Outer diameter: 45 mm

Part no.	Unit
8500220	1 roll/25 m



IMP

Thermal paper LPQ80, 80 m

For GSI TympStar™ Pro.
Width: 80 mm
Inner diameter: 13 mm
Outer diameter: 80 mm

Part no.	Unit
8509893	1 roll/80 m



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 rolls

ABR IMP OAE

API300 thermal paper for Titan™, 10 m

Width: 57 mm
Inner diameter: 12 mm
Outer diameter: 35 mm

Part no.	Unit
8500230	1 roll/10 m



ABR IMP OAE

Sanibel™ MPT-III thermal paper, 20 m

For Sanibel™ MPT-III printer which is supplied with Interacoustics® AT235™.

Width: 79 mm
Inner diameter: 12 mm
Outer diameter: 45 mm

Part no.	Unit
8500610	1 roll/20 m



IMP

Adhesive thermal paper LPQ80, 195 labels

For GSI TympStar™ Pro.
Width: 82 mm
Inner diameter: 25 mm
Outer diameter: 88 mm

Part no.	Unit
8509894	1 roll/195 labels



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 rolls

IMP OAE

Thermal paper for Sanibel™ MPT-II printer, 20 m

For the Sanibel™ MPT-II printer which is supplied with MAICO® ERO-SCAN™, Grason-Stadler® Corti™ and Interacoustics® OtoRead™.

Width: 57 mm
Inner diameter: 9 mm
Outer diameter: 40 mm

Part no.	Unit
8029305	1 roll/20 m



AUD IMP

TPR26™ thermal paper for Interacoustics®, 25 m

Fits most popular audiometers and tympanometers.

Width: 112 mm
Inner diameter: 12 mm
Outer diameter: 45 mm

Part no.	Unit
8500210	1 roll/25 m



Heine/
KaWe
100 pcs

Welch
Allyn
250 pcs

No sharp
edges



AUD IMP

Otoscope specula

Disposable ear speculums for Heine/KaWe and Welch Allyn otoscopes size 2.5 and 4 mm.

Part no.	Description	Unit
8507618	4 mm Speculum Economy for Heine/KaWe	100 pcs
8507594	2.5 mm Speculum Economy for Heine/KaWe	100 pcs
8514367	4 mm Speculum for Heine/KaWe	100 pcs
8514368	2.5 mm Speculum for Heine/KaWe	100 pcs
8514369	4 mm Speculum Genuine for Heine/KaWe	50 pcs
8514370	2.5 mm Speculum Genuine for Heine/KaWe	50 pcs
8107750	2.5 mm Speculum for Welch Allyn	250 pcs
8108252	4 mm Speculum for Welch Allyn	250 pcs



AUD IMP

HEINE® Mini3000® Otoscope

Xenon Halogen Technology. Swiveling viewing window. Minimal reflection. On/off switch in clip.

Part no.	Unit
8108833	1 pc



No sharp
edges

50 pcs

AUD BAL IMP

Disposable speculas for video otoscope

Fits the AURICAL™ OTOcam 300, Interacoustics® VIOT and XION™ video otoscopes.

Part no.	Unit
8500470	Bag/50 pcs

2

LATEX



No sharp
edges

100 pcs

AUD BAL IMP

Video Oscopes specula

Fits the new Interacoustics® Viot and OriVision otoscope.

Part no.	Description	Unit
8515083	KaWe 2.5 mm speculum	Bag/100 pcs
8511728	KaWe 4 mm speculum	Bag/100 pcs

2

LATEX



5 pcs

AUD IMP

EARNimals®

Cartoon faces that are placed on the otoscope. Make ear examinations less scary for kids.

Part no.	Unit
8011501	Bag/5 pcs (assorted)



HIF

Ear coupler sealing wax

Use to create a seal between coupler and hearing instrument.


Part no.	Unit
8500080	1 pc



ABR AUD OAE

Alkaline battery 1.5V AA

Part no.	Unit
8011226	1 pc



HIF

Battery adaptor

For Interacoustics Affinity™, Callisto™, Equinox™, MedRx Avant Hit+ and Siemens UNITY 2.

Part no.	Description	Unit
8011198	BAA675	1 pc
8011199	BAA13	1 pc
8011200	BAA312	1 pc
8011201	BAA10	1 pc
8011202	BAA5	1 pc



IMP HIF

Lithium battery 3.7V

For Interacoustics® Callisto™ and Titan™.

Part no.	Unit
8009978	1 pc



HIF

Dreve® mold syringe – clear plastic

Easy to fill, inject and clean. Provides a stable and easy injection.

Part no.	Unit
8102613	1 pc





Sanibel^{*} Supply

Please order at your local special
instrument distributor:



Insert Foam Tips 3A
Adult/Medium



Sanibel™ Infant EarCup™
with snap electrodes



Sanibel™ ADI ear tips,
8 mm mushroom



Insert Foam Tips 3B
Pediatric/Small



Sanibel™ REM probe tubes,
type A



Ambu® Blue Sensor™
N-00-S

Disclaimer

This catalogue has been published for worldwide distribution. Some parts described may NOT BE AVAILABLE in your area. Your local distributor will be pleased to give you details. For advice on compatibility and instructions of use, please consult your operation manual or local distributor. Sanibel Supply A/S assumes no legal liability or responsibility for the accuracy, typing errors or mistakes, correctness or completeness of any information in this flyer. The information in this flyer was correct to the best of our knowledge at the time of printing. Sanibel Supply A/S reserves the right to make changes at any time in specification of products described in this flyer without notice or obligation. MADSEN AccuScreen® and MADSEN OTOflex® are registered trademarks of Otometrics A/S. Echo-Screen®, AuDx®, Navigator® Pro, and ABAer® are registered trademarks of Natus Medical Incorporated. Sanibel™ is a registered trademark of Sanibel Supply®. All other Instrument brands or model names are copyrighted or trademarked by their respective owners and are used only for reference.

