

# ENT microscopes

Karl Kaps and Zeiss





Since 1983, **Heinemann Medizintechnik** is one of the leading equipment providers in **ENT medical technology**.

The company is not only a **manufacturer** of high-quality ENT treatment units, but also a **competent** partner for practice planning as well as your sales **partner** for complete ENT equipment in **hospitals** and private practices.

# ENT microscopes



**SOM 22**

**SOM 62**

## **SOM model range**

Standard equipment:

- 3-level magnification changer
- Lens  $f = 250$  mm
- Binocular straight tube  $f = 125$  mm
- 2x Eyepieces for spectacle wearers, WF 12.5xV or WF 10xV
- LED illumination unit

# Examples of additional equipment



Microscope holder with fine focussing



Rapid lens changer



Beam splitter with endoscope adapter



Lenses



Handles



45° coupler



Lens fine focus



Monitor holder



Vario objective



Observer tubes

Ordering Information	REF
KAPS SOM 22 (installation in unit) with coaxial <b>LED</b> cold light illumination	535-0500005
KAPS SOM 62 (floor installation) with cold light illumination	535-0500014
Beam splitter with endoscope adapter	535-101890000
Lens f = 200 mm	535-101630000
Lens f = 250 mm	535-101640000
Lens f = 300 mm	535-101660000
Lens f = 400 mm	535-101670000
Lens fine focus with integrated lens f = 200 mm, focus area 15 mm	535-101701000
Vario lens f = 200-400 mm	535-162700000
Cold light adapter with green filter	535-162335000
Rapid lens changer for 3x changer (without lens)	535-102181000
Magnification 5x step changer	535-101120200
Rapid lens changer for 5x changer (without lens)	535-102191000
Wall bracket for wall mounting	535-102110000
T-handle, infinitely adjustable	535-102185100
Monitorhalterung am SOM-62	535-162205000

Additional equipment available on request

# ZEISS OPMI pico

Comprising the OPMI pico with 5x changer, integrated green and orange filter, a T-handle and integrated cold light fibreoptic illumination with rapid lamp changer

- Straight wide-field tube 45° f = 170 mm
- 2 eyepieces with 10x magnification
- Lens f = 250 mm
- Wall, ceiling or floor stand including illumination unit and fibreoptic light cable
- 2 sleeves for tubes with screw thread
- Dust cover
- 6-pack of asepsis caps 12 mm





Ordering information	REF
Zeiss Opmi pico (installation on unit)	547-0500074
Endoport (beam splitter with endoscope adapter)	547-1153971

# Contact

G. Heinemann Medizintechnik GmbH

Leibnizstraße 13-15

24568 Kaltenkirchen

Germany

Phone: +49 (0)4191 - 95379-0

Fax: +49 (0)4191 - 9537937

Email: [info@heinemann-ent.de](mailto:info@heinemann-ent.de)

Web: [www.heinemann-ent.de](http://www.heinemann-ent.de)

